

# Unitek™ Gemini Bracket

## MBT Versatile+ Appliance System

MAXILLARY	TORQUE	ANG	IN/OUT		L/R	APC™ II Adhesive Coated Appliance System	
			in.	mm		.018 in.	.022 in.
			Central	17°		4°	0.03
	17°	4°			R	3 119-614	3 119-714
Lateral	10°	8°	0.04	1.06	L	3 119-615	3 119-715
	10°	8°			R	3 119-616	3 119-716
Cuspid**	-7°	8°	0.03	0.80	L	3 119-671	3 119-771
	-7°	8°			R	3 119-672	3 119-772
Cuspid w/hook**	-7°	8°	0.03	0.80	L	3 119-673	3 119-773
	-7°	8°			R	3 119-674	3 119-774
Cuspid*	0°	8°	0.03	0.80	L	3 119-675	3 119-775
	0°	8°			R	3 119-676	3 119-776
Cuspid w/hook*	0°	8°	0.03	0.80	L	3 119-677	3 119-777
	0°	8°			R	3 119-678	3 119-778
1st/2nd Bicuspid	-7°	0°	0.03	0.83	L/R	3 119-635	3 119-735
1st/2nd Bicuspid w/hook	-7°	0°	0.03	0.83	L	3 119-665	3 119-765
	-7°	0°			R	3 119-666	3 119-766
2nd Bicuspid	-7°	0°	0.04	1.06	L/R	3 119-679	3 119-779
2nd Bicuspid w/hook	-7°	0°	0.04	1.06	L	3 119-681	3 119-781
	-7°	0°			R	3 119-680	3 119-780

MANDIBULAR	TORQUE	ANG	IN/OUT		L/R	APC™ II Adhesive Coated Appliance System	
			in.	mm		.018 in.	.022 in.
			Anterior	-6°		0°	0.05
Cuspid**	-6°	3°	0.03	0.83	L	3 119-683	3 119-783
	-6°	3°			R	3 119-684	3 119-784
Cuspid w/hook**	-6°	3°	0.03	0.83	L	3 119-685	3 119-785
	-6°	3°			R	3 119-686	3 119-786
Cuspid*	0°	3°	0.03	0.83	L	3 119-687	3 119-787
	0°	3°			R	3 119-688	3 119-788
Cuspid w/hook*	0°	3°	0.03	0.83	L	3 119-689	3 119-789
	0°	3°			R	3 119-690	3 119-790
1st Bicuspid	-12°	2°	0.03	0.83	L	3 119-691	3 119-791
	-12°	2°			R	3 119-692	3 119-792
1st Bicuspid w/hook	-12°	2°	0.03	0.83	L	3 119-693	3 119-793
	-12°	2°			R	3 119-694	3 119-794
2nd Bicuspid	-17°	2°	0.03	0.83	L	3 119-695	3 119-795
	-17°	2°			R	3 119-696	3 119-796
2nd Bicuspid w/hook	-17°	2°	0.03	0.83	L	3 119-697	3 119-797
	-17°	2°			R	3 119-698	3 119-798

\*\*Denotes which cuspid brackets are included in the kits.

3 Use this prefix to order APC™ II Adhesive Coated Appliance System

### Kits

	APC™ II Adhesive Coated Appliance System	
	.018 in.	.022 in.
One Case Kit 3 with hook*	3 119-143	3 119-142
One Case Kit 3, 4, 5 with hook*	3 119-145	3 119-144
One Case Kit 3 with hook**	3 119-153	3 119-152
One Case Kit 3, 4, 5 with hook**	3 119-155	3 119-154

\*\*Denotes which cuspid brackets are included in the kits.

3 Use this prefix to order APC™ II Adhesive Coated Appliance System

**3M Unitek**  
**EMEA Headquarter**  
 Ohmstr. 3  
 86899 Landsberg / Germany  
 www.3MUnitek.net

For additional information please  
 contact your local representative.

Please recycle.  
 © 2010 3M. All rights reserved.  
 016-929 EMEA 1005

**3M Unitek**