

Unitek™ Gemini Brackets

Roth Prescription*

MAXILLARY	TORQUE	ANG	IN/OUT		L/R	APC™ II Adhesive Coated Appliance System	
			in.	mm		.018 in.	.022 in.
Central	14°	5°	0.031	0.79	L	3 119-621	3 119-721
	14°	5°			R	3 119-622	3 119-722
Lateral	7°	8°	0.040	1.02	L	3 119-623	3 119-723
	7°	8°			R	3 119-624	3 119-724
Cuspid	0°	10°	0.029	0.74	L	3 119-625	3 119-725
	0°	10°			R	3 119-626	3 119-726
Cuspid w/hook	0°	10°	0.029	0.74	L	3 119-661	3 119-761
	0°	10°			R	3 119-662	3 119-762
1st/2nd Bicuspid	-7°	0°	0.032	0.81	L/R	3 119-635	3 119-735
1st/2nd Bicuspid w/hook	-7°	0°	0.032	0.81	L	3 119-665	3 119-765
	-7°	0°			R	3 119-666	3 119-766

MANDIBULAR	TORQUE	ANG	IN/OUT		L/R	APC™ II Adhesive Coated Appliance System	
			in.	mm		.018 in.	.022 in.
Anterior	-1°	0°	0.049	1.24	L/R	3 119-640	3 119-740
Cuspid	-11°	6°	0.032	0.81	L	3 119-611	3 119-711
	-11°	6°			R	3 119-612	3 119-712
Cuspid w/hook	-11°	6°	0.032	0.81	L	3 119-659	3 119-759
	-11°	6°			R	3 119-660	3 119-760
Cuspid	0°	6°	0.030	0.76	L	3 119-645	3 119-745
	0°	6°			R	3 119-646	3 119-746
Cuspid w/hook	0°	6°	0.030	0.76	L	3 119-663	3 119-763
	0°	6°			R	3 119-664	3 119-764
1st Bicuspid	-17°	0°	0.032	0.81	L/R	3 119-655	3 119-755
1st Bicuspid w/hook	-17°	0°	0.032	0.81	L	3 119-667	3 119-767
	-17°	0°			R	3 119-668	3 119-768
2nd Bicuspid	-22°	0°	0.032	0.81	L/R	3 119-656	3 119-756
2nd Bicuspid w/hook	-22°	0°	0.032	0.81	L	3 119-669	3 119-769
	-22°	0°			R	3 119-670	3 119-770

*3M Unitek version of this prescription. No endorsement by the Doctor is implied.

3 Use this prefix to order APC™ II Adhesive Coated Appliance System

Kits

	APC™ II Adhesive Coated Appliance System	
	.018 in.	.022 in.
One Case Kit 3 with hook	3 119-162	3 119-172
One Case Kit 3 with hook**	3 119-182	3 119-192
One Case Kit 3, 4, 5 with hook	3 119-163	3 119-173
One Case Kit 3, 4, 5 with hook**	3 119-183	3 119-193

**The kits marked with two asterisks include the -11° torque lower cuspid.

3 Use this prefix to order APC™ II Adhesive Coated Appliance System

3M Unitek
EMEA Headquarter
 Ohmstr. 3
 86899 Landsberg / Germany
 www.3MUnitek.net

For additional information please
 contact your local representative.

Please recycle.
 © 2010 3M. All rights reserved.
 016-932 EMEA 1005

3M Unitek